



## SW1 Complete Weight loss and Skin Toning Machine



### Specifications

**Vacuum liposuction probe (40 KHZ cavitation with vacuum suction)**

Used for larger areas of stubborn fat.

**Ultra sound probe (40 KHZ cavitation)**

Used for softer areas of fat.

**Ultra therapy probes x 3 sizes (1MHZ cavitation)**

Used for smaller areas of fat and cellulite - e.g. face, arms etc.

**Waveform Bi Polar Probe (RF)**

Used for all over skin tightening and toning.



## Operation

### Main Selection Screen




### RF- Bipolar RF Interface



**Bipolar treatment heads.**

Bi-polar treatment head: 15mm, 32mm  
 RF output Levels: 8  
 RF energy Level: 0~100  
 Frequency: 1MHz / 10MHz / 1MHz & 10MHz

	Measurement	Main parts to be used for
1#	φ15	Smaller part, such as the eye, the neck, etc.
2#	φ32	Bigger part, such as the face, etc.

**Bipolar parameters**

The below parameters are only for Bipolar RF machine. You should select the right parameters for different clients according to their body state, age and the skin conditions.

Treatment head	Treatment Area	RF output Mode	RF energy Level	Instructions
Φ15	Eyes and small areas	± 2	±20J	Output mode level 3 should not exceed 3 for sensitive areas.
		± 5	±15J	Suitable for normal skin.
Φ32	Medium sized parts such as face, neck and arms	± 3	±20J	Suitable for sensitive areas
		± 6	±25J	Suitable for normal skin

**Application of RF**

1. Facial skin tightening
2. Fine wrinkle removal



3. The mandible and the neck skin tightening
4. Improve the conditions of sagging skin
5. Forehead wrinkle removal
6. Breast lifting
7. Renew back contours
8. Renew belly contours
9. Renew buttocks and legs contours
10. Hand care
11. Improving the channel between the lip and the nose
12. Weight losing of the arms and legs

## Fat Cavitation



## Main technical parameter

1. Ultra Cavitation: Frequency 40K Hz or 1MHz
2. Energy density: 0~5W/cm<sup>2</sup>
3. Output energy: 100W(40K) or 30W(1M)
4. Size of hand piece: 40K:  $\Phi$ 50mm and 1M:  $\Phi$ 8mm,  $\Phi$ 45mm,  $\Phi$ 65mm

## 40KHZ Cavitation Interface



Body  Logy®



**Liposuction 40K plus Suction**



**1M Cavitation  
1M Interface**

